

检验报告书 / Certificate of Analysis

产品名称 Product Name	HEK293 CDM15	检验编号 Analysis ID	20220712001A
产品编号 Catalog No.	DP15	报告日期 Report Date	2022.07.28
生产批号 Lot Number	DP15220712001A	贮存条件 Storage Requirement	2 ~ 8°C 避光 Protected from light
生产日期 MFG Date	2022.07.12	有效日期 Expiration Date	2024.07.11
检验依据 Test Standard	检验方法标准操作规程 Sop for inspection methods		

检验项目 Test Categories	标准要求 Test Specifications	检验结果 Test Results
外观 Physical Appearance	白色至淡粉色及相近颜色粉末 White to pale pink and similar color power	白色至淡粉色及相近颜色粉末 White to pale pink and similar color power
澄清度 Clarity	澄清 Clarify	合格 Pass
pH 值 pH Value	6.8 ~ 7.4	合格 Pass
渗透压 Osmolality	270 ~ 320 (mOsmol/Kg)	合格 Pass
细菌内毒素 Endotoxin	< 10 EU/mL	合格 Pass
微生物限度 Microbial Limit	需氧菌 Aerobic bacteria: <200CFU/g	合格 Pass
	霉菌酵母菌 Molds and yeasts: <50CFU/g	合格 Pass

结论: 本品按成品检验操作规程检测, 结果均已记录。

Conclusion: This product has been tested according to the finished product inspection operating procedures and the results have been recorded.

QC 复核人: 黄晶晶
Quality Control:

QA 复核人: [Signature]
Quality Assurance:

审核人: [Signature]
Approver:

