

## 检验报告书 / Certificate of Analysis

产品名称 Product Name	PDE RecTrypsin-LQ03(10×)	检验编号 Analysis ID	230620PDELQ03001B
产品编号 Catalog No.	PDE-LQ03	报告日期 Report Date	2023.07.04
生产批号 Lot Number	PDE-LQ03230620001.B	贮存条件 Storage Requirement	2 ~ 8°C 避光 Protected from light
生产日期 MFG Date	2023.06.20	有效日期 Expiration Date	2024.06.19
检验依据 Test Standard	检验方法标准操作规程 Sop for inspection methods		

检验项目 Test Categories	标准要求 Test Specifications	检验结果 Test Results
外观 Physical Appearance	无色澄清液体 Colorless transparent liquid	无色澄清液体 Colorless transparent liquid
澄清度 Clarity	澄清 Clear	澄清 Clear
pH 值 pH Value	6.9 ~ 7.4	合格 PASS
渗透压 Osmolality	270 ~ 320 (mOsmol/Kg)	合格 PASS
酶活力 Activity Assay (BAEE)	≥1600U/mL	2534.7 U/mL
体外生物活性 In-Vitro Bioassay	符合 Acceptable	合格 PASS
细菌内毒素 Endotoxin	<1.0 EU/mL	合格 PASS
无菌检查 Sterility check	阴性 Negative	合格 PASS

结论: 本品按成品检验操作规程检测, 结果均已记录。

Conclusion: This product has been tested according to the finished product inspection operating procedures and the results have been recorded.

QC 复核人: 陈义梅  
Quality Control:

QA 复核人: 付紫娟  
Quality Assurance:

