

产品规格 / PRODUCT SPECIFICATIONS
---HEK293 CDM15

产品名称 Product Name	HEK293 CDM15 With Glucose (6.00g/L) Without L-Glutamine
产品编号 Product Code	DP15
产品特点 Characteristics	<ul style="list-style-type: none"> ● 化学成分限定，无动物源性成分，适用于HEK293、HEK293F、293T等悬浮高密度培养及腺病毒及AAV(Adeno-Associated Viral Vector)腺相关病毒转染表达，适合大规模生产应用。 ● Chemical defined media, animal origin-free medium specially designed for HEK293, HEK293F, 293T cell line. It can be used for the growth of HEK293 suspension cells high-density culture adenovirus and AAV transfection and expression, suitable for large-scale manufacturing applications.
形式 Format	粉末，非无菌 / Powder, Non-sterile
包装规格 Pack Size	50L/100L (20.30 g/L)
存放条件 Storage	2°C ~ 8°C，避光 / 2°C ~ 8°C, Protected from light
有效期 Shelf life	24 个月 / 24 Months