

**产品规格 / PRODUCT SPECIFICATIONS**
**---CHO CDM29**

|                                       |                                                                                                                                                                                                                                                                                                                                                         |
|---------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>产品名称</b><br><b>Product Name</b>    | CHO CDM29<br>With Glucose (6.00g/L)<br>Without L-Glutamine                                                                                                                                                                                                                                                                                              |
| <b>产品编号</b><br><b>Product Code</b>    | DP29                                                                                                                                                                                                                                                                                                                                                    |
| <b>产品特点</b><br><b>Characteristics</b> | <ul style="list-style-type: none"> <li>● 化学成分限定，无动物源组分，无蛋白，适合CHO系列细胞悬浮高密度大规模培养，可提高蛋白及抗体的表达与质量。</li> <li>● Chemical defined, animal origin-free, protein-free medium optimized for the growth of high-density and large-scale suspension cultures of CHO series cells, and can improve the expression and quality of antibodies and proteins.</li> </ul> |
| <b>形式</b><br><b>Format</b>            | 粉末，非无菌 / Powder, Non-sterile                                                                                                                                                                                                                                                                                                                            |
| <b>包装规格</b><br><b>Pack Size</b>       | 50L/100L (21.12 g/L)                                                                                                                                                                                                                                                                                                                                    |
| <b>存放条件</b><br><b>Storage</b>         | 2°C ~ 8°C，避光 / 2°C ~ 8°C, Protected from light                                                                                                                                                                                                                                                                                                          |
| <b>有效期</b><br><b>Shelf life</b>       | 24 个月 / 24 Months                                                                                                                                                                                                                                                                                                                                       |