

产品规格 / PRODUCT SPECIFICATIONS
---MDCK SFM01

产品名称 Product Name	MDCK SFM01 With Glucose (8.20g/L), L-Glutamine
产品编号 Product Code	DP17
产品特点 Characteristics	<ul style="list-style-type: none"> ● 无血清、无动物源，适用于MDCK细胞悬浮高密度大规模培养,适应于疫苗高效表达。 ● Serum-free, animal origin-free medium suitable for the high-density culture of MDCK, suitable for vaccine large-scale manufacturing applications.
形式 Format	粉末，非无菌 / Powder, Non-sterile
包装规格 Pack Size	50L/100L (27.21 g/L)
存放条件 Storage	2°C ~ 8°C，避光 / 2°C ~ 8°C, Protected from light
有效期 Shelf life	24 个月 / 24 Months