

产品规格 / **PRODUCT SPECIFICATIONS**

--- **DMEM (Low Glucose)**

<p>产品名称 Product Name</p>	<p>DMEM (Low Glucose) With Glucose (1.00g/L), L-Glutamine, Phenol red, Sodium pyruvate Without NaHCO₃, HEPES</p>
<p>产品编号 Product Code</p>	<p>DP23</p>
<p>产品特点 Characteristics</p>	<ul style="list-style-type: none"> ● 支持多种哺乳动物细胞系的生长。 ● Basal medium for supporting the growth of many different mammalian cell lines.
<p>形式 Format</p>	<p>粉末，非无菌 / Powder, Non-sterile</p>
<p>包装规格 Pack Size</p>	<p>50L/100L (9.99 g/L)</p>
<p>存放条件 Storage</p>	<p>2°C ~ 8°C, 避光 / 2°C ~ 8°C, Protected from light</p>
<p>有效期 Shelf life</p>	<p>36 个月 / 36 Months</p>