

产品规格 / PRODUCT SPECIFICATIONS

---细胞保存液 (MCP CryoSpr-DHQ 10(A)) Cryoprotectant (MCP CryoSpr-DHQ 10(A))

产品名称 Product Name	细胞保存液 (MCP CryoSpr-DHQ 10(A)) Cryoprotectant (MCP CryoSpr-DHQ 10(A))
产品编号 Product Code	CP004-100
产品特点 Characteristics	<ul style="list-style-type: none">● 用于多种不同类型细胞的保存及运送, 如外周血单个核细胞、造血干细胞、间充质干细胞、免疫细胞(T、NK、DC等) 以及CHO、HEK293、293T、JurKat及K562等。● It is used for the preservation and transport of many different types of cells, such as peripheral blood mononuclear cells, hematopoietic stem cells, mesenchymal stem cells, immune cells (T, NK, DC, etc.) and CHO, HEK293, 293T, JurKat and K562 etc.
形式 Format	液体, 无菌 / Liquid, Asepsis
包装规格 Pack Size	100mL/瓶
存放条件 Storage	2°C ~ 8°C, 避光 / 2°C ~ 8°C, Protected from light
有效期 Shelf Life	24个月/24 Months