

**产品规格 / PRODUCT SPECIFICATIONS**
**--- DPBS**

<b>产品名称</b> <b>Product Name</b>	DPBS Without Ca <sup>2+</sup> , Mg <sup>2+</sup>
<b>产品编号</b> <b>Product Code</b>	DP12
<b>产品特点</b> <b>Characteristics</b>	<ul style="list-style-type: none"> <li>● Dulbecco 磷酸盐缓冲溶液 (DPBS) 是一种平衡盐溶液, 适用于各种细胞培养应用。</li> <li>● Dulbecco's phosphate-buffered saline (DPBS) is a balanced salt solution used for a variety of cell culture applications.</li> </ul>
<b>形式</b> <b>Format</b>	粉末, 非无菌 / Powder, Non-sterile
<b>包装规格</b> <b>Pack Size</b>	50L/100L (9.55 g/L)
<b>存放条件</b> <b>Storage</b>	2°C ~ 8°C, 避光 / 2°C ~ 8°C, Protected from light
<b>有效期</b> <b>Shelf life</b>	36 个月 / 36 Months