

产品规格 / **PRODUCT SPECIFICATIONS**

--- 重组胰蛋白酶消化液
Recombinant Trypsin Solution

产品名称 Product Name	重组胰蛋白酶消化液 PDE RecTrypsin-LQ02(1×)
产品编号 Product Code	PDE-LQ02
产品特点 Characteristics	<ul style="list-style-type: none"> ● 不含动物源性，可替代猪或牛胰蛋白酶，用于免疫细胞治疗、疫苗、药物筛选、抗体等领域的细胞消化过程。本产品含EDTA，不含酚红。 ● animal origin-free, can replace porcine or bovine trypsin, and is used in cell digestion processes in immunocell therapy, vaccines, drug screening, antibodies and other fields. This product contains EDTA and does not contain phenol red.
形式 Format	液体，无菌 / Liquid, sterile
包装规格 Pack Size	100, 500mL/瓶
存放条件 Storage	2°C ~ 8°C, 避光 / 2°C ~ 8°C, Protected from light
有效期 Shelf life	12个月 / 12 Months