

产品规格 / PRODUCT SPECIFICATIONS
**--- 全能核酸酶
BenzoNuclease**

产品名称 Product Name	全能核酸酶 BenzoNuclease
产品编号 Product Code	210615TO, 5, 25, 50, 5000KU
产品特点 Characteristics	<ul style="list-style-type: none"> ● 可用于在科学研究中降低细胞上清和细胞裂解液的粘度，提高蛋白纯化效率及功能研究，也可以应用在基因治疗、病毒纯化、疫苗生产、蛋白和多糖类制药工业作为宿主残留核酸去除试剂。 ● It can be used in scientific research to reduce the viscosity of cell supernatant and cell lysate, improve the efficiency and function of protein purification, and can also be used in gene therapy, virus purification, vaccine production, protein and polysaccharide pharmaceutical industry as host residual nucleic acid removal reagent.
形式 Format	液体，无菌 / Liquid, sterile
包装规格 Pack Size	5KU/20μL, 50KU/200μL, 500KU/2mL, 5000KU/20mL
存放条件 Storage	- 25 ~ - 15°C
有效期 Shelf life	24个月 / 24 Months