

产品规格 / PRODUCT SPECIFICATIONS

---NK细胞无血清培养基试剂盒 NK Serum-Free Medium Kit

产品名称 Product Name	NK细胞无血清培养基试剂盒 NK Serum-Free Medium Kit		
产品编号 Product Code	CM-NKKIT03		
产品特点 Characteristics	<ul style="list-style-type: none">● 无血清，无动物源组分，适合NK高效快速培养。● Serum-Free, animal origin-free, optimized for the NK high efficiency and rapid expansion.		
形式 Format	液体，无菌 / Liquid, Asepsis		
包装规格 Pack Size	2.2L		
存放条件 Storage	NK细胞无血清扩增培养基 NK serum-free expansion medium	CM022A	2 ~ 8°C
	NK细胞无血清激活培养基 NK serum-free activating medium	CM022B	2 ~ 8°C
	包被因子	NK-CF	-20°C
	激活因子	NK-AF	-20°C
	扩增因子	NK-EF	-20°C
	血清处理剂	NK-STR	2 ~ 8°C
有效期 Shelf Life	12 个月/12 Months		