

产品规格 / **PRODUCT SPECIFICATIONS**

---Tobitec FedVigor02

产品名称 Product Name	Tobitec FedVigor02 With Glucose (75.00g/L) Without L-Glutamine
产品编号 Product Code	DP47
产品特点 Characteristics	● 化学成分确定的补料培养基，无动物源性成分，无蛋白。可应用于CHO及HEK293的瞬时转染表达，有效提高细胞密度及产物表达量。 ● Chemical defined, animal origin-free, protein-free medium. Tobitec FedVigor01 can be applied to the transient transfection expression of CHO and HEK293 cells, effectively improving the cell density and product expression.
形式 Format	粉末，非无菌 / Powder, Non-sterile
包装规格 Pack Size	5L/10L（195.41 g/L）
存放条件 Storage	2°C ~ 8°C，避光 / 2°C ~ 8°C, Protected from light
有效期 Shelf life	24 个月 / 24 Months