

产品规格 / PRODUCT SPECIFICATIONS

---HEK Trans kit01

产品名称 Product Name	HEK293瞬转试剂盒 HEK293 Transient Transfection Kit
产品编号 Product Code	HEK Trans kit01
产品特点 Characteristics	● HEK293瞬转试剂盒化学成分限定，适用于不同亚型的HEK293细胞高密度培养及瞬时转染表达，各组分配合使用可显著提高转染性能和蛋白表达水平，具有优良的稳定性和重现性。 ● The HEK293 Transient Kit is chemically component-limited and is suitable for high-density culture and transient transfection expression of different subtypes of HEK293 cells. The combined use of each group can significantly improve the transfection performance and protein expression level, and has excellent stability and reproducibility.
形式 Format	液体，无菌 / Liquid, sterile
包装规格 Pack Size	HEK293 CDM26 Pro(26-2)- 1000ml×1 Tobitec FedVigor02 - 125ml×2 CHE STE01 - 20ml×1 Tobitec TR01 - 2ml×2
存放条件 Storage	2°C ~ 8°C，避光 / 2°C ~ 8°C, Protected from light
有效期 Shelf life	HEK293 CDM26 Pro(26-2): 12个月 / 12 Months Tobitec FedVigor02: 6个月 / 6 Months CHE STE01: 6个月 / 6Months Tobitec TR01: 6个月 / 6 Months