Version 2，2026

产品规格 / PRODUCT SPECIFICATIONS

---HEK293 CDM26 Pro

产品名称 **Product Name**	HEK293 CDM26 Pro With Glucose (7.00g/L)，L-Glutamine
产品编号 **Product Code**	DP26-2
产品特点 **Characteristics**	● 化学成分限定,适应于HEK293（293F、293T、293H、Expi293等）悬浮高密度培养，蛋白及慢病毒等转染表达，适合大规模生产应用。 ● Chemical defined media, animal origin-free medium specially designed for HEK293.It can be used for the growth of HEK293（293F,293T,293H,Expi293） suspension cells high-density culture, lentivirus and proteins expression, suitable for large-scale manufacturing applications
形式 **Format**	粉末，非无菌 / Powder, Non-sterile
包装规格 **Pack Size**	50L/100L（23.00 g/L）
存放条件 **Storage**	2℃ ~ 8℃，避光 / 2℃ ~ 8℃, Protected from light
有效期 **Shelf life**	24 个月 / 24 Months

Shanghai Tofflon Bio-reagent Co.,Ltd No.139 Duhui Road,Minhang District,Shanghai,China http://bitec. tofflon.com Tel: +86 21-64909996-393 Fax: +86 21-64909996-730