Version 3，2026

产品规格 / PRODUCT SPECIFICATIONS

---CHO CDM10

产品名称 **Product Name**	CHO CDM10 With Glucose (7.00g/L), L-Glutamine
产品编号 **Product Code**	DP10
产品特点 **Characteristics**	● 化学成分限定，无动物源组分，无蛋白，适合CHO悬浮高密度大规模培养，尤其是CHO-S及CHO-K1，可提高蛋白及抗体的表达与质量。 ● Chemical defined, animal origin-free, protein-free medium optimized for the growth of high-density and large-scale suspension cultures of CHO (especially for CHO-S and CHO-K1) cells, and can improve the expression and quality of antibodies and proteins.
形式 **Format**	粉末，非无菌 / Powder, Non-sterile
包装规格 **Pack Size**	50L/100L（23.00 g/L）
存放条件 **Storage**	2℃ ~ 8℃，避光 / 2℃ ~ 8℃, Protected from light
有效期 **Shelf life**	24 个月 / 24 Months

Shanghai Tofflon Bio-reagent Co.,Ltd No.139 Duhui Road,Minhang District,Shanghai,China http://bitec. tofflon.com Tel: +86 21-64909996-393 Fax: +86 21-64909996-730