产品规格 / PRODUCT SPECIFICATIONS

---MDCK SFM01

产品名称 **Product Name**	MDCK SFM01 With Glucose (8.20g/L), L-Glutamine
产品编号 **Product Code**	DP17
产品特点 **Characteristics**	● 无血清、无动物源，适用于MDCK细胞悬浮高密度大规模培养,适应于疫苗高效表达。 ● Serum-free, animal origin-free medium suitable for the high-density culture of MDCK, suitable for vaccine large-scale manufacturing applications.
形式 **Format**	粉末，非无菌 / Powder, Non-sterile
包装规格 **Pack Size**	50L/100L（27.21 g/L）
存放条件 **Storage**	2℃ ~ 8℃，避光 / 2℃ ~ 8℃, Protected from light
有效期 **Shelf life**	24 个月 / 24 Months