产品规格 / PRODUCT SPECIFICATIONS

--- BHK SFM01

产品名称 **Product Name**	BHK SFM01 With Glucose (8.70g/L), L-Glutamine
产品编号 **Product Code**	DP16
产品特点 **Characteristics**	● 无动物源、无血清培养基，适合BHK21细胞无血清生长培养。 ● Animal origin-free, serum-free medium suitable for the culture of BHK21 cell.
形式 **Format**	粉末，非无菌 / Powder, Non-sterile
包装规格 **Pack Size**	50L/100L（24.724 g/L）
存放条件 **Storage**	2℃ ~ 8℃，避光 / 2℃ ~ 8℃, Protected from light
有效期 **Shelf life**	24 个月 / 24 Months