产品规格 / PRODUCT SPECIFICATIONS

--- LMH SFM01

产品名称 **Product Name**	LMH SFM01 With Glucose (7.20g/L), L-Glutamine, Sodium pyruvate Without $NaHCO_3$, Phenol red, HEPES
产品编号 **Product Code**	DP42
产品特点 **Characteristics**	● 无动物源、无血清培养基，适合 LMH 细胞无血清生长培养。 ● Animal origin-free, serum-free medium suitable for the culture of LMH1 cell. ● 支持 LMH 细胞高效扩增禽腺病毒。 ● Support high-efficiency amplification of avian adenovirus.
形式 **Format**	粉末，非无菌 / Powder, Non-sterile
包装规格 **Pack Size**	50L/100L（23.30 g/L）
存放条件 **Storage**	2℃ ~ 8℃，避光 / 2℃ ~ 8℃, Protected from light
有效期 **Shelf life**	24 个月 / 24 Months