产品规格 / PRODUCT SPECIFICATIONS

--- Leibovitz's L-15

产品名称 **Product Name**	Leibovitz's L-15 With L-Galactose，L-Glutamine， Phenol red，Sodium pyruvate Without Glucose，$NaHCO_3$，HEPES
产品编号 **Product Code**	DP07
产品特点 **Characteristics**	● 可支持 HEP-2 猴肾脏细胞以及胚胎和成人组织的原代组织块细胞的生长。 ● It can support the growth of HEP-2 monkey kidney cells as well as primary tissue block cells in embryonic and adult tissues.
形式 **Format**	粉末，非无菌 / Powder, Non-sterile
包装规格 **Pack Size**	50L/100L（14.08 g/L）
存放条件 **Storage**	2℃ ~ 8℃，避光 / 2℃ ~ 8℃, Protected from light
有效期 **Shelf life**	36 个月 / 36 Months