Tofflon Bio-reagent

Tobitec®

# 产品规格 / PRODUCT SPECIFICATIONS

## --- MEM

| | |
|---|---|
| **产品名称**<br>**Product Name** | MEM<br>With Glucose (1.00g/L), L-Glutamine, Phenol red<br>Without $NaHCO_3$, Sodium pyruvate, HEPES |
| **产品编号**<br>**Product Code** | DP01 |
| **产品特点**<br>**Characteristics** | ● 支持多种哺乳动物细胞系的生长。<br>● Basal medium for supporting the growth of many different mammalian cell lines. |
| **形式**<br>**Format** | 粉末，非无菌 / Powder, Non-sterile |
| **包装规格**<br>**Pack Size** | 50L/100L（9.53 g/L） |
| **存放条件**<br>**Storage** | 2℃ ~ 8℃，避光 / 2℃ ~ 8℃, Protected from light |
| **有效期**<br>**Shelf life** | 36 个月 / 36 Months |