产品规格 / PRODUCT SPECIFICATIONS

---Insect SFM01

产品名称 **Product Name**	Insect SFM01 With Glucose (11.00g/L), L-Glutamine
产品编号 **Product Code**	DP11
产品特点 **Characteristics**	● 无血清,无动物源组分，专门为Insect细胞系开发，可用于Insect悬浮高密度培养即杆状病毒转染表达，适用于大规模生产应用。 ● 联合 Insect TE01使用。 ● No serum, no animal origin components, specially developed for Insect cell line, can be used for suspension high density culture of Insect namely baculovirus transfection expression, suitable for large-scale production application. ● Using combined with Insect TE01.
形式 **Format**	粉末，非无菌 / Powder, Non-sterile
包装规格 **Pack Size**	100L（41.60 g/L）
存放条件 **Storage**	2℃ ~ 8℃，避光 / 2℃ ~ 8℃, Protected from light
有效期 **Shelf life**	36 个月 / 36 Months