产品规格 / PRODUCT SPECIFICATIONS

--- RPMI 1640

产品名称 **Product Name**	RPMI 1640 With Glucose (2.00g/L), L-Glutamine, Phenol red Without $NaHCO_3$, Sodium pyruvate, HEPES
产品编号 **Product Code**	DP03
产品特点 **Characteristics**	● 支持多种哺乳动物细胞系的生长。 ● Basal medium for supporting the growth of many different mammalian cell lines.
形式 **Format**	粉末，非无菌 / Powder, Non-sterile
包装规格 **Pack Size**	50L/100L（10.41 g/L）
存放条件 **Storage**	2℃ ~ 8℃，避光 / 2℃ ~ 8℃, Protected from light
有效期 **Shelf life**	36 个月 / 36 Months