产品规格 / PRODUCT SPECIFICATIONS

---Insect TE01

产品名称 **Product Name**	Insect TE01
产品编号 **Product Code**	TE01
产品特点 **Characteristics**	● 无血清，无动物源组分。 ● Insect SFM01 添加剂（Insect TE01经乙醇溶解后按1mL/L添加）。 ● Serum-free, animal origin-free. ● Additive Insect SFM01 (Add 1mL /L after Insect TE01 dissolved in ethanol)
形式 **Format**	悬浊液，非无菌/Suspension, Non-sterile
包装规格 **Pack Size**	3.2g/瓶（适配 100 L Insect SFM01 培养基）/ 3.2 g/bottle (Compatible with 100 L Insect SFM01 medium)
存放条件 **Storage**	-20℃，避光 / -20℃, Protected from light
有效期 **Shelf life**	18 个月 / 18 Months