Version 2，2026

产品规格 / PRODUCT SPECIFICATIONS

---Ham's F-12K

产品名称 **Product Name**	Ham's F-12K With Glucose (1.26g/L), L-Glutamine, Phenol red, Sodium pyruvate Without $NaHCO_3$， HEPES
产品编号 **Product Code**	DP38
产品特点 **Characteristics**	● 可支持原代人类肝细胞以及减血清环境中的部分大鼠和鸡肝细胞生长。 ● It can support the growth of primary human hepatocytes and some rat and chicken hepatocytes in a reduced serum environment.
形式 **Format**	粉末，非无菌 / Powder, Non-sterile
包装规格 **Pack Size**	50L/100L（11.09 g/L）
存放条件 **Storage**	2℃ ~ 8℃，避光 / 2℃ ~ 8℃, Protected from light
有效期 **Shelf life**	36 个月 / 36 Months

Shanghai Tofflon Bio-reagent Co.,Ltd No.139 Duhui Road,Minhang District,Shanghai,China http://bitec. tofflon.com Tel: +86 21-64909996-393 Fax: +86 21-64909996-730