Version 2，2026

产品规格 / PRODUCT SPECIFICATIONS

---Tobitec FedVigor01

产品名称 **Product Name**	Tobitec FedVigor01 With Glucose (75.00g/L) Without L-Glutamine
产品编号 **Product Code**	DP46
产品特点 **Characteristics**	● 化学成分确定的补料培养基，无动物源性成分，无蛋白。可应用于CHO及HEK293的稳定转染表达，有效提高细胞密度及产物表达量。 ● Chemical defined, animal origin-free, protein-free medium. Tobitec FedVigor01 can be applied to the stable transfection expression of CHO and HEK293 cells, effectively improving the cell density and product expression.
形式 **Format**	粉末，非无菌 / Powder, Non-sterile
包装规格 **Pack Size**	5L/10L（196.23 g/L）
存放条件 **Storage**	2℃ ~ 8℃，避光 / 2℃ ~ 8℃, Protected from light
有效期 **Shelf life**	24 个月 / 24 Months

Shanghai Tofflon Bio-reagent Co.,Ltd No.139 Duhui Road,Minhang District,Shanghai,China http://bitec. tofflon.com Tel: +86 21-64909996-393 Fax: +86 21-64909996-730