Version 1，2026

产品规格 / PRODUCT SPECIFICATIONS

---Ecoli CDM02

产品名称 **Product Name**	Ecoli CDM02 With Glucose (16.5g/L)
产品编号 **Product Code**	DPM02
产品特点 **Characteristics**	● 化学成分限定，无动物源组分，无蛋白，适合大肠杆菌的高密度悬浮培养，可实现优质质粒的高效生产。 ● Chemical defined, animal origin-free, protein-free medium optimized for the growth of high-density of *E. coli*, and can achieve efficient production of high-quality plasmids.
形式 **Format**	粉末，非无菌 / Powder, Non-sterile
包装规格 **Pack Size**	50L/100L（35.5g/L）
存放条件 **Storage**	2℃ ~ 8℃，避光 / 2℃ ~ 8℃, Protected from light
有效期 **Shelf life**	18 个月 / 18 Months

Shanghai Tofflon Bio-reagent Co.,Ltd No.139 Duhui Road,Minhang District,Shanghai,China http://bitec. tofflon.com Tel: +86 21-64909996-393 Fax: +86 21-64909996-730